



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Fluorescent Dye 750-I Streptavidin ABN-U0304

Article Name	Fluorescent Dye 750-I Streptavidin
Biozol Catalog Number	ABN-U0304
Supplier Catalog Number	U0304
Alternative Catalog Number	ABN-U0304-200
Manufacturer	Abnova
Category	Sonstiges
Application	Conjugation, FC
Product Description	Fluorescence Dye 750-I with Streptavidin conjugate. Compared with Alexa Fluor 750, Fluorescent Dye 750-I is much brighter with absorption at 633nm. Fluorescence emission of Fluorescence Dye 750-I is well separated from that of other commonly used red
Buffer	Water
Application Dilute	Labeling via Streptavidin for fluorescence imaging and other fluorescence-based biochemical analysis. For most fluorescent streptavidin applications, a concentration of 1-20 ug/mL is sufficient. The optimal working dilution should be determined by the