

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

HotTaq DNA Polymerase ABN-D0011

Article Name	HotTaq DNA Polymerase
Biozol Catalog Number	ABN-D0011
Supplier Catalog Number	D0011
Alternative Catalog Number	ABN-D0011-500,ABN-D0011-1000,ABN-D0011-2500
Manufacturer	Abnova
Category	Molekularbiologie
Application	PCR
Product Description	HotTaq DNA polymerase catalyzes the 5-3 synthesis of DNA. The enzyme has been proven not to have the 3 - 5 exonuclease activity. Before the first PCR step, the HotTaq DNA polymerase should be activated by 15 minute incubation at 95-97C....
Concentration	10 units/uL
Buffer	20 mM Tris-HCl (pH 8.0), 1 mM DTT, 0.1 mM EDTA, 100 mM KCl, 0.5% Nonidet P40, 0.5% Tween 20 and 50% glycerol.