



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Load 6x Loading Dye Bromophenol Blue ABN-D0030

Article Name	Load 6x Loading Dye Bromophenol Blue
Biozol Catalog Number	ABN-D0030
Supplier Catalog Number	D0030
Alternative Catalog Number	ABN-D0030-6
Manufacturer	Abnova
Category	Sonstiges
Application	EP
Product Description	The 6x loading dye contains bromophenol blue and xylene cyanol FF for easy visual tracking of DNA migration during electrophoresis. The ready-to-use solution can be used for loading samples on agarose or polyacrylamide gels
Application Dilute	1/6 volume of the Load 6x Loading Dye solution should be added to the DNA sample.