



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## HepG2/C3A (human hepatoblastoma) nuclear extract lysate (nondenatured), Human ABN-L019V4

Article Name	HepG2/C3A (human hepatoblastoma) nuclear extract lysate (non-denatured), Human
Biozol Catalog Number	ABN-L019V4
Supplier Catalog Number	L019V4
Alternative Catalog Number	ABN-L019V4-50
Manufacturer	Abnova
Host	Human
Category	Zellen/Zellkultur
Application	IP, WB
Species Reactivity	Human
Product Description	Nuclear extract cell lysate (non-denatured)
Buffer	In Buffer D.
Source	Liver
Application Dilute	Use it directly for immuno-precipitation, or heat lysate with SDS gel loading buffer to 95C for 5 minutes followed by rapid cooling for western blot application. If dissociating conditions are required, add reducing agent prior to heating.