



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rat kidney tissue lysate (non-denatured) ABN-L039W2

Article Name	Rat kidney tissue lysate (non-denatured)
Biozol Catalog Number	ABN-L039W2
Supplier Catalog Number	L039W2
Alternative Catalog Number	ABN-L039W2-200
Manufacturer	Abnova
Host	Rat
Category	Zellen/Zellkultur
Application	IP, WB
Species Reactivity	Rat
Product Description	Whole tissue lysate (non-denatured)
Buffer	In modified RIPA Lysis Buffer.
Source	Kidney
Application Dilute	Use it directly for immuno-precipitation, or heat lysate with SDS gel loading buffer to 95C for 5 minutes followed by rapid cooling for western blot application. If dissociating conditions are required, add reducing agent prior to heating.