



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

HEK293/Human NAPI-IIb Stable Cell Line ABS-CHEK-ATP116-2VIALS

Article Name	HEK293/Human NAPI-IIb Stable Cell Line
Biozol Catalog Number	ABS-CHEK-ATP116-2VIALS
Supplier Catalog Number	CHEK-ATP116-2Vials
Alternative Catalog Number	ABS-CHEK-ATP116-2VIALS
Manufacturer	AcroBiosystems
Category	Zellen/Zellkultur
Species Reactivity	Human
Product Description	NAPIIIb, NaPi2b, Npt2b, or SLC34A2 enables sodium:phosphate symporter activity. Acts upstream of or within in utero embryonic development and phosphate ion transport. Located in apical plasma membrane and brush border. Is expressed in several structu
Source	NAPI-IIb, NaPi-2b, Npt2b, or SLC34A2 enables sodium:phosphate symporter activity. Acts upstream of or within in utero embryonic development and phosphate ion transport. Located in apical plasma membrane and brush border. Is expressed in several structure
Expression System	HEK293
Form	Liquid
Target	NAPI-IIb
Application Notes	Useful for cell-based NAPI-IIb binding assay