

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## ActiveMax Human CD19 muBeads, premium grade (for cells) ABS-MBS-C002-2.5MG

Article Name	ActiveMax Human CD19 muBeads, premium grade (for cells)
Biozol Catalog Number	ABS-MBS-C002-2.5MG
Supplier Catalog Number	MBS-C002-2.5mg
Alternative Catalog Number	ABS-MBS-C002-2.5MG
Manufacturer	AcroBiosystems
Host	Human
Category	Zellen/Zellkultur
Species Reactivity	Human
Product Description	ActiveMax Human CD19 $\mu$ Beads, premium grade (for cells) are uniform, superparamagnetic beads of 5.5 $\mu$ m with streptavidin on its surface and coupled with biotinylated Human CD19 Protein, expressed from human 293 cells (HEK293) and contains AA Pro 20 Ly
Source	ActiveMax Human CD19 $\mu$ Beads, premium grade (for cells) are uniform, superparamagnetic beads of 5.5 $\mu$ m with streptavidin on its surface and coupled with biotinylated Human CD19 Protein, expressed from human 293 cells (HEK293) and contains AA Pro 20 - Lys
Form	Powder
Target	CD19

## Application Notes

ActiveMax Human CD19  $\mu$ Beads, premium grade (for cells) are designed to stimulate in vitro CD19-specific CAR-T cells or UCAR-T cells, similar to the tumor cell lines that express human CD19 antigen. It can be used as follows: Evaluating the characteristi