

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

ActiveMax Human MSLN muBeads, premium grade (for cells) ABS-MBS-C003-2.5MG

Article Name	ActiveMax Human MSLN muBeads, premium grade (for cells)
Biozol Catalog Number	ABS-MBS-C003-2.5MG
Supplier Catalog Number	MBS-C003-2.5mg
Alternative Catalog Number	ABS-MBS-C003-2.5MG
Manufacturer	AcroBiosystems
Host	Human
Category	Zellen/Zellkultur
Species Reactivity	Human
Product Description	ActiveMax Human MSLN μ Beads, premium grade (for cells) are uniform, superparamagnetic beads of 5.5 μ m with streptavidin on its surface and coupled with biotinylated Human MSLN Protein, expressed from human 293 cells (HEK293) and contains AA Glu 296 G
Source	ActiveMax Human MSLN μ Beads, premium grade (for cells) are uniform, superparamagnetic beads of 5.5 μ m with streptavidin on its surface and coupled with biotinylated Human MSLN Protein, expressed from human 293 cells (HEK293) and contains AA Glu 296 - Gl
Form	Powder
Target	Mesothelin

Apı	plication	Notes
, VP	piication	110103

ActiveMax Human MSLN μ Beads, premium grade (for cells) are designed to stimulate in vitro MSLN-specific CAR-T cells or UCAR-T cells, similar to the tumor cell lines that express human MSLN antigen. It can be used as follows: Evaluating the characteristi