

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## ActiveMax Human Her2 muBeads, premium grade (for cells) ABS-MBS-C006-10MG

Article Name	ActiveMax Human Her2 muBeads, premium grade (for cells)
Biozol Catalog Number	ABS-MBS-C006-10MG
Supplier Catalog Number	MBS-C006-10mg
Alternative Catalog Number	ABS-MBS-C006-10MG
Manufacturer	AcroBiosystems
Host	Human
Category	Zellen/Zellkultur
Species Reactivity	Human
Product Description	ActiveMax Human Her2 $\mu$ Beads, premium grade (for cells) are uniform, superparamagnetic beads of 5.5 $\mu$ m with streptavidin on its surface and coupled with biotinylated Human Her2 Protein, expressed from human 293 cells (HEK293) and contains AA Thr 23 Th
Source	ActiveMax Human Her2 $\mu$ Beads, premium grade (for cells) are uniform, superparamagnetic beads of 5.5 $\mu$ m with streptavidin on its surface and coupled with biotinylated Human Her2 Protein, expressed from human 293 cells (HEK293) and contains AA Thr 23 - Thr
Form	Powder
Target	Her2

## **Application Notes**

ActiveMax Human Her2  $\mu$ Beads, premium grade (for cells) are designed to stimulate in vitro Her2-specific CAR-T cells or UCAR-T cells, similar to the tumor cell lines that express human Her2 antigen. It can be used as follows: Evaluating the characteristi