

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

ActiveMax Streptavidin muBeads, premium grade (for cells) ABS-MBS-C009-100MG

Article Name	ActiveMax Streptavidin muBeads, premium grade (for cells)
Biozol Catalog Number	ABS-MBS-C009-100MG
Supplier Catalog Number	MBS-C009-100mg
Alternative Catalog Number	ABS-MBS-C009-100MG
Manufacturer	AcroBiosystems
Category	Zellen/Zellkultur
Product Description	ActiveMax Streptavidin coupledµbeads, premium grade (for cells) are produced under sterile manufacturing conditions (ISO 5), and no animal or humanderived components are used throughout the production process. It is produced under our rigorous qualit
Source	ActiveMax Streptavidin coupled-µbeads, premium grade (for cells) are produced under sterile manufacturing conditions (ISO 5), and no animal- or human-derived components are used throughout the production process. It is produced under our rigorous qualit
Form	Powder
Target	Streptavidin Beads
Application Notes	ActiveMax Streptavidin coupled-µbeads, premium grade (for cells)can be used for coupling of any biotinylated proteins or antibodies, and are ideal for numerous applications as follows: When coupled with biotinylated CAR target proteins, it can be used t