

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## ActiveMax Human DLL3 muBeads, premium grade (for cells) ABS-MBS-C018-10MG

Article Name	ActiveMax Human DLL3 muBeads, premium grade (for cells)	
Biozol Catalog Number	ABS-MBS-C018-10MG	
Supplier Catalog Number	MBS-C018-10mg	
Alternative Catalog Number	ABS-MBS-C018-10MG	
Manufacturer	AcroBiosystems	
Host	Human	
Category	Zellen/Zellkultur	
Species Reactivity	Human	
Product Description	ActiveMax Human DLL3 $\mu$ Beads, premium grade (for cells) are uniform, superparamagnetic beads of 5.5 $\mu$ m with streptavidin on its surface and coupled with biotinylated Human DLL3 Protein, expressed from human 293 cells (HEK293) and contains AA Ala 27 Le	
Source	ActiveMax Human DLL3 $\mu$ Beads, premium grade (for cells) are uniform, superparamagnetic beads of 5.5 $\mu$ m with streptavidin on its surface and coupled with biotinylated Human DLL4 Protein, expressed from human 293 cells (HEK293) and contains AA Ala 27 - Leu	
Form	Powder	
Target	DLL3	

Αp	plicatio	n Notes
, vp	piioatio	11110100

ActiveMax Human DLL3  $\mu$ Beads, premium grade (for cells) are designed to stimulate in vitro DLL3-specific CAR-T cells or UCAR-T cells, similar to the tumor cell lines that express human DLL3 antigen. It can be used as follows: Evaluating the characteristi