

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Anti-SARS-CoV-2 Spike RBD Antibody-coupled Magnetic Beads ABS-MBS-K014-2MG

Article Name	Human Anti-SARS-CoV-2 Spike RBD Antibody-coupled Magnetic Beads
Biozol Catalog Number	ABS-MBS-K014-2MG
Supplier Catalog Number	MBS-K014-2mg
Alternative Catalog Number	ABS-MBS-K014-2MG
Manufacturer	AcroBiosystems
Category	Antikörper
Product Description	The Human AntiSARSCoV2 Spike RBD Antibodycoupled Magnetic Beads is produced by coupling biotinylated antiRBD antibody to streptavidinconjugated magnetic beads. The precoupled beads are ready to use for capturing spike protein from your sample with hi...
Source	The Human Anti-SARS-CoV-2 Spike RBD Antibody-coupled Magnetic Beads is produced by coupling biotinylated anti-RBD antibody to streptavidin-conjugated magnetic beads. The pre-coupled beads are ready to use for capturing spike protein from your sample with
Form	Powder
Target	Spike RBD
Application Notes	This product is intended for immunocapture, biopanning and flow cytometry. This product is produced non-sterile.