

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

CF770 Goat Anti-Rabbit IgG (H+L), highly cross-absorbed, 2 mg/mL BOT-20078-1

Article Name	CF770 Goat Anti-Rabbit IgG (H+L), highly cross-absorbed, 2 mg/mL
Biozol Catalog Number	BOT-20078-1
Supplier Catalog Number	20078-1
Alternative Catalog Number	BOT-20078-1-50UL
Manufacturer	Biotium
Category	Antikörper
Conjugation	CF770
Product Description	This product is prepared by labeling highly cross-adsorbed goat anti- rabbit IgG (H+L) with CF770 dye. CF dyes are Biotiums line of next- generation fluorescent dyes with advantages in brightness, photostability, and conjugate specificity compared to o
Concentration	2 mg/mL
Buffer	PBS, 50% glycerol, 2 mg/ml BSA, 0.05% azide
Source	Animal