

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## CF640R Goat Anti-Rabbit IgG (H+L), Highly Cross-Adsorbed, 2 mg/mL BOT-20176-1

Article Name	CF640R Goat Anti-Rabbit IgG (H+L), Highly Cross-Adsorbed, 2 mg/mL
Biozol Catalog Number	BOT-20176-1
Supplier Catalog Number	20176-1
Alternative Catalog Number	BOT-20176-1-50UL
Manufacturer	Biotium
Category	Antikörper
Conjugation	CF640R
Product Description	This product is prepared by labeling highly cross-adsorbed goat anti- rabbit IgG (H+L) with CF640R dye. CF dyes are Biotiums line of next- generation fluorescent dyes with advantages in brightness, photostability, and conjugate specificity compared to
Concentration	2 mg/mL
Buffer	PBS, 50% glycerol, 2 mg/ml BSA, 0.05% azide
Source	Animal