

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## CF770 Donkey Anti-Guinea Pig IgG (H+L), Highly Cross-Adsorbed, 2 mg/mL BOT-20242

Article Name	CF770 Donkey Anti-Guinea Pig IgG (H+L), Highly Cross-Adsorbed, 2 mg/mL
Biozol Catalog Number	BOT-20242
Supplier Catalog Number	20242
Alternative Catalog Number	BOT-20242-250UL
Manufacturer	Biotium
Category	Antikörper
Conjugation	CF770
Product Description	This product is prepared by labeling highly cross-adsorbed donkey anti-guinea pig IgG (H+L) with CF770 dye. CF dyes are Biotiums line of next-generation fluorescent dyes with advantages in brightness, photostability, and conjugate specificity compare
Concentration	2 mg/mL
Buffer	PBS, 50% glycerol, 2 mg/ml BSA, 0.05% azide
Source	Animal