

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

CF660C Goat Anti-Rabbit IgG (H+L), Highly Cross-Adsorbed, 2 mg/mL BOT-20369-500UL

Article Name	CF660C Goat Anti-Rabbit IgG (H+L), Highly Cross-Adsorbed, 2 mg/mL
Biozol Catalog Number	BOT-20369-500UL
Supplier Catalog Number	20369-500uL
Alternative Catalog Number	BOT-20369-500UL-500UL
Manufacturer	Biotium
Category	Antikörper
Conjugation	CF660C
Product Description	This product is prepared by labeling highly cross-adsorbed goat anti- rabbit IgG (H+L) with CF660C dye. CF dyes are Biotiums line of next- generation fluorescent dyes with advantages in brightness, photostability, and conjugate specificity compared to
Concentration	2 mg/mL
Buffer	PBS, 50% glycerol, 2 mg/ml BSA, 0.05% azide
Source	Animal