

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## CF532 Goat Anti-Rabbit IgG (H+L), Highly Cross-Adsorbed BOT-20469-1MG

Article Name	CF532 Goat Anti-Rabbit IgG (H+L), Highly Cross-Adsorbed
Biozol Catalog Number	BOT-20469-1MG
Supplier Catalog Number	20469-1mg
Alternative Catalog Number	BOT-20469-1MG-1MG
Manufacturer	Biotium
Category	Antikörper
Conjugation	CF532
Product Description	This product is prepared by labeling highly cross-adsorbed goat anti- rabbit IgG (H+L) with CF532 dye. CF dyes are Biotiums line of next- generation fluorescent dyes with advantages in brightness, photostability, and conjugate specificity compared to o
Concentration	N/A (lyophilized)
Buffer	PBS, 15 mg/mL BSA, 20 mg/mL trehalose (after reconstitution)
Source	Animal