



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SAP-Exo Kit BYT-ORB1733626

Article Name	SAP-Exo Kit
Biozol Catalog Number	BYT-ORB1733626
Supplier Catalog Number	orb1733626
Alternative Catalog Number	BYT-ORB1733626-500
Manufacturer	Biorbyt
Category	Sonstiges
Product Description	The Kit contains two hydrolytic enzymes, recombinant Shrimp Alkaline Phosphatase (rSAP) and Exonuclease I (Exo I). The combination of these enzymes ensures complete dephosphorylation of dNTPs and degradation of residual primers, enabling downstream
Application Notes	SAP-Exo Kit is a quick, easy and reliable enzymatic cleanup reagent for PCR product cleanup. Removes excess primers and dNTPs, Fast 20 minute protocol, Add directly to PCR product, 100% Sample Recovery, Scalable for different reaction sizes, No interfere