

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

ERG (Vascular Endothelial & Prostate Marker) Antibody, Clone: [rERG/6843], Unconjugated, Recombinant BYT-ORB2313971

| | |
|----------------------------|---|
| Article Name | ERG (Vascular Endothelial & Prostate Marker) Antibody, Clone: [rERG/6843], Unconjugated, Recombinant |
| Biozol Catalog Number | BYT-ORB2313971 |
| Supplier Catalog Number | orb2313971 |
| Alternative Catalog Number | BYT-ORB2313971-100 |
| Manufacturer | Biorbyt |
| Category | Antikörper |
| Species Reactivity | Human |
| Immunogen | Recombinant fragment (around aa450 to 479) (C-terminus) of human ERG (exact sequence is proprietary) |
| Conjugation | Unconjugated |
| Product Description | ERG (Vascular Endothelial & Prostate Marker) Antibody... |
| Clonality | Recombinant |
| Concentration | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml |
| Clone Designation | [rERG/6843] |
| Molecular Weight | 55kDa |
| UniProt | P11308 |

| | |
|--------------------|---|
| Application Dilute | Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT)(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95 C followed by cooling at RT for 20 minutes)Optimal dilution fo |
| Application Notes | Positive Control: Human tonsil, Ewing sarcoma or prostate adenocarcinoma tissues. |