



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Seleno-D,L-ethionine USB-S0065-10

Article Name	Seleno-D,L-ethionine
Biozol Catalog Number	USB-S0065-10
Supplier Catalog Number	S0065-10
Alternative Catalog Number	USB-S0065-10-50
Manufacturer	US Biological
Category	Biochemikalien
Product Description	Used as a substrate for methionine adenosyltransferase from rat liver. Melting Point: 247-250C Solubility: Hydrochloric Acid (aq), Boiling Methanol, Hot Water, Storage and Stability: May be stored at RT for short-term only. Long-term storage is recom
Form	Cream Colored Shiny Crystals