

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## SAMS Peptide (Substrate for AMP-activated Protein Kinase) USB-S0090

Article Name	SAMS Peptide (Substrate for AMP-activated Protein Kinase)
Biozol Catalog Number	USB-S0090
Supplier Catalog Number	S0090
Alternative Catalog Number	USB-S0090-250
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	The AMP-activated protein kinase (AMPK) cascades appear to be sensors for cellular energy status. The members of the AMP-activated/SNF-1-related subfamily of protein kinases are thought to be pivotal components of these highly conserved protein kinas
Molecular Weight	1.8
NCBI	000429
UniProt	Q00266
Purity	95%
Form	Supplied as a liquid in sterile ddH2O.