



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SAP-I (Speract) USB-S0096-01

Article Name	SAP-I (Speract)
Biozol Catalog Number	USB-S0096-01
Supplier Catalog Number	S0096-01
Alternative Catalog Number	USB-S0096-01-5
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Sequence (linear): Gly-Phe-Asp-Leu-Asn-Gly-Gly-Gly-Val-Gly Storage and Stability: Lyophilized powder may be stored at 4C for short-term only. Stable for 12 months at -20C. Reconstitute to nominal volume (see reconstitution instructions for peptides)
Purity	~95%. Purified by HPLC.
Form	Supplied as a lyophilized powder.