



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Streptavidin (TRITC) USB-S7973-81

Article Name	Streptavidin (TRITC)
Biozol Catalog Number	USB-S7973-81
Supplier Catalog Number	S7973-81
Alternative Catalog Number	USB-S7973-81-1
Manufacturer	US Biological
Category	Molekularbiologie
Application	ELISA, IF, WB
Product Description	Label: Rhodamine (Tetramethylrhodamine isothiocyanate) (MW: 444D) Absorption Wavelength: 550 nm Emission Wavelength: 570 nm Fluorochrome/Protein Ratio: 2.7 moles TRITC per mole of Streptavidin Applications: Suitable for use in ELISA (Fluorescent), We
Purity	Purified by affinity chromatography. IEP Assay shows a single precipitin arc against anti-steptavidin. No reaction was observed against anti-avidin.
Form	Supplied as a lyophilized powder from PBS, pH 7.2, 10mg/ml BSA (IgG and Protease free), 0.01% sodium azide. Labeled with TRITC. Reconsitute with 1ml sterile ddH2O.