



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Streptavidin (FITC) USB-S7973-96A

Article Name	Streptavidin (FITC)
Biozol Catalog Number	USB-S7973-96A
Supplier Catalog Number	S7973-96A
Alternative Catalog Number	USB-S7973-96A-1
Manufacturer	US Biological
Category	Molekularbiologie
Application	FLISA, IF, WB
Product Description	Streptavidin/FITC can be used in conjunction with biotinylated antibodies as an avidin/biotin labeling system for flow cytometry. Streptavidin (MW $\sim$ 66kD) binds specifically with biotin (244D). It is derived from the bacterium Streptomyces avidinii an
Purity	Purified by ion exchange chromatography. IEP Assay shows a single precipitin arc against source serum protein and IgG.
Form	Supplied as a lyophilized powder in PBS, pH 7.2, 10mg/ml BSA, 0.09% sodium azide.