

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Tissue Plasminogen Activator (tPA), ATA-FPR-CMK Inactivated, Human USB-T5600-07S

Article Name	Tissue Plasminogen Activator (tPA), ATA-FPR-CMK Inactivated, Human
Biozol Catalog Number	USB-T5600-07S
Supplier Catalog Number	T5600-07S
Alternative Catalog Number	USB-T5600-07S-100
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Prepared from human tPA by incubation with an excess of ATA- Phe-Pro-Arg-CMK. After labeling at the active site, remaining free peptide is removed by dialysis. Source: Human sc-tPA, ATA-FPR-CMK Inactivated Recommended Dilution: Optimal dilutions to b
Purity	95% (SDS-PAGE)
Form	Supplied as a liquid in 0.4M Hepes, 0.1M sodium chloride, pH 7.4.