



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anhydrotrypsin, Bovine USB-T8656-20

Article Name	Anhydrotrypsin, Bovine
Biozol Catalog Number	USB-T8656-20
Supplier Catalog Number	T8656-20
Alternative Catalog Number	USB-T8656-20-100
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Treatment of trypsin with phenylmethylsulfonyl fluoride (PMSF) followed by a strong base leads to the conversion of Serine 195 to dehydroalanine at the active site (1). Alpha and beta forms of anhydrotrypsin are separated by immobilized soybean tryps
Molecular Weight	83
Purity	~95% (SDS-PAGE)
Form	Supplied as a liquid in 1mM HCl, 10mM CaCl2.