

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## ZO-2 (Tight Junction-associated Polypeptide, Zonula Occludens 2) (MaxLight 650), Rabbit USB-Z0760-ML650

Article Name	ZO-2 (Tight Junction-associated Polypeptide, Zonula Occludens 2) (MaxLight 650), Rabbit
Biozol Catalog Number	USB-Z0760-ML650
Supplier Catalog Number	Z0760-ML650
Alternative Catalog Number	USB-Z0760-ML650-100
Manufacturer	US Biological
Host	Rabbit
Category	Antikörper
Application	FLISA, IHC, WB
Immunogen	Synthetic peptide corresponding to the central portion of the canine ZO-2 protein.
Product Description	MaxLight(TM) 650 is a new Far-IR stable dye conjugate comparable to Alexa Fluor(TM)647, DyLight(TM)649, Cy5(TM) and offers better labeling efficiency, brighter imaging and increased immunodetection. Absorbance (655nm), Emission (676nm), Extinction Co
Purity	Purified by immunoaffinity chromatography.
Form	Supplied as a liquid in PBS, pH 7.2. No preservatives added. Labeled with MaxLight(TM)650.