

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rabbit Anti Dog IgG (H&L)-AbFluor 680 AST-SA1013

Article Name	Rabbit Anti Dog IgG (H&L)-AbFluor 680
Biozol Catalog Number	AST-SA1013
Supplier Catalog Number	SA1013
Alternative Catalog Number	AST-SA1013-100, AST-SA1013-1
Manufacturer	Assay Biotechnology
Host	Rabbit
Category	Antikörper
Application	FC, IF, WB
Species Reactivity	Canine
Conjugation	AbFluor 680
Product Description	Assay Biotechnology secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offe
Source	Polyclonal, Rabbit IgG
Purity	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Form	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Di	lute
----------------	------

Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR Odyssey equipment would be recommended.