

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Human Neurotrophin-3 (NTF3) BYT-ORB1096095

Recombinant Human Neurotrophin-3 (NTF3)
BYT-ORB1096095
orb1096095
BYT-ORB1096095-20, BYT-ORB1096095-100, BYT-ORB1096095-1
Biorbyt
Proteine/Peptide
Recombinant Human Neurotrophin-3(NTF3)
17.5 kDa
P20783
If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Homo sapiens (Human)
Greater than 85% as determined by SDS-PAGE.
Liquid or Lyophilized powder
YAEHKSHRGEYSVCDSESLWVTDKSSAIDIRGHQVTVLGEIKTGNSPVKQYFY ETRCKEARPVKNGCRGIDDKHWNSQCKTSQTYVRALTSENNKLVGWRWIRI DTSCVCALSRKIGRT

Application	Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0~mg/mL.We recommend to add 5-50% of glycerol (final concentration