

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Phleum pratense Pollen allergen Phl p 6 (PHLPVI) BYT-ORB1096116

Article Name	Recombinant Phleum pratense Pollen allergen Phl p 6 (PHLPVI)
Biozol Catalog Number	BYT-ORB1096116
Supplier Catalog Number	orb1096116
Alternative Catalog Number	BYT-ORB1096116-20, BYT-ORB1096116-100, BYT-ORB1096116-1
Manufacturer	Biorbyt
Category	Proteine/Peptide
Product Description	Recombinant Phleum pratense Pollen allergen Phl p 6(PHLPVI)...
Molecular Weight	19.2 kDa
UniProt	<a href="#">P43215</a>
Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Source	Phleum pratense (Common timothy)
Purity	Greater than 90% as determined by SDS-PAGE.
Form	Liquid or Lyophilized powder
Sequence	GKATTEEQKLIEDVNASFRAAMATTANVPPADKYKTFEAAFTVSSKRNLADAV SKAPQLVPKLDEVYNAAYNAADHAAPEDKYEAFVLHFSEALRIIAGTPEVHAVK PGA

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration