

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Bovine coronavirus Non-structural protein 2a (2a) BYT-ORB1096642

Article Name	Recombinant Bovine coronavirus Non-structural protein 2a (2a)
Biozol Catalog Number	BYT-ORB1096642
Supplier Catalog Number	orb1096642
Alternative Catalog Number	BYT-ORB1096642-20, BYT-ORB1096642-100, BYT-ORB1096642-1
Manufacturer	Biorbyt
Category	Proteine/Peptide
Product Description	Recombinant Bovine coronavirus Non-structural protein 2a(2a)...
Molecular Weight	33.2 kDa
UniProt	Q80B88
Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Source	Bovine coronavirus (strain Mebus) (BCoV) (BCV)
Purity	Greater than 90% as determined by SDS-PAGE.
Form	Liquid or Lyophilized powder
Sequence	MAVAYADKPNHFINFPLTQFEGFVLNYKGLQFQLLDEGVDCIKQTAPHISLAML DIQPEDYRSVDVAIQEVIDDMHWGEGFQIKFDNPLILGRCIVLDVKGVEELHD DLVNYIRDKGCVADQSRKWIGHCTIAQLTDAALSIKENVDFINSMQFNKITIN PSSPARLEIVKLGAEEKDGFYETIVSHWMGIRFEYNPPTDKLAMIMGYCCLEVV RKELEEGDLPENDDDAWFKLSYHYENNSWFFRHVYRKSS

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration