

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Bovine coronavirus Non-structural protein 2a (2a) BYT-ORB1096658

| | |
|----------------------------|--|
| Article Name | Recombinant Bovine coronavirus Non-structural protein 2a (2a) |
| Biozol Catalog Number | BYT-ORB1096658 |
| Supplier Catalog Number | orb1096658 |
| Alternative Catalog Number | BYT-ORB1096658-20, BYT-ORB1096658-100, BYT-ORB1096658-1 |
| Manufacturer | Biorbyt |
| Category | Proteine/Peptide |
| Product Description | Recombinant Bovine coronavirus Non-structural protein 2a(2a)... |
| Molecular Weight | 34.4 kDa |
| UniProt | Q8V438 |
| Buffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Source | Bovine coronavirus (strain 98TXSF-110-LUN) (BCoV-LUN) (BCV) |
| Purity | Greater than 85% as determined by SDS-PAGE. |
| Form | Liquid or Lyophilized powder |
| Sequence | MAVAYADKPNHFINFPLTQFQGFVLNYKGLQFQLLDEGVDCCKIQTAPHISLAM LDIQPEDYRSVDVAIQEVIDDMHWGEGFQIKFENPHILGRCIVLDVKGVEELH DDLVNYYIRDKGCVDDQSRKWIGHCTIAQLTDAALSIGKGVDFINSMQFNKITI NPSSPARLEIVKLGAEKKDGFYETIASHWMGIRFEYNPPTDKLAMIMGYCCLEV VRKELEEGDLPENDDAWFKLSYHYENNSWFFRHVYRKSS |

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration