

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Sheep anti Rabbit Glyceraldehyde-3-phosphate dehydrogenase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE066/BIO

Article Name	Sheep anti Rabbit Glyceraldehyde-3-phosphate dehydrogenase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE066/BIO
Supplier Catalog Number	NE066/Bio
Alternative Catalog Number	NMB-NE066/BIO
Manufacturer	NordicMubio
Host	Sheep
Category	Antikörper
Conjugation	Biotin
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and
Clonality	Monoclonal
Clone Designation	[Polyclonal]
NCBI	100009074
UniProt	P46406
Buffer	Biotin-coupled hyperimmune sheep IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile dis

Source	Glyceraldehyde-3-phosphate dehydrogenase isolated and purified from rabbit muscle. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Biotin-coupled hyperimmune sheep IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting.