

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Sheep anti Rabbit Glycerol-3-phosphate dehydrogenase, Clone: [Polyclonal], Monoclonal NMB-NE070/PAB

Article Name	Sheep anti Rabbit Glycerol-3-phosphate dehydrogenase, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE070/PAB
Supplier Catalog Number	NE070/PAb
Alternative Catalog Number	NMB-NE070/PAB
Manufacturer	NordicMubio
Host	Sheep
Category	Antikörper
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and
Clonality	Monoclonal
Clone Designation	[Polyclonal]
UniProt	P08507
Buffer	Purified hyperimmune sheep IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitute

## BIOZOL

Source	Glycerol-3-phosphate dehydrogenase isolated and purified from rabbit muscle. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Purified hyperimmune sheep IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	Immunoprecipitation.