

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Sheep anti Mouse IgG Fab, Clone: [Polyclonal], Monoclonal NMB-SHAM/FAB

Article Name	Sheep anti Mouse IgG Fab, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-SHAM/FAB
Supplier Catalog Number	ShAM/Fab
Alternative Catalog Number	NMB-SHAM/FAB
Manufacturer	NordicMubio
Host	Sheep
Category	Antikörper
Product Description	The antiserum is reacting with the Fab subunit of intact IgG, IgA and IgM and other Ig classes of both light chain types, with their Fab or F(ab)2 subunits and with free light chains of kappa and lambda type as tested in immunoelectrophoresis and dou
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Delipidated, heat inactivated lyophilized stable whole serum No preservative added. Total protein and IgG concentration in the antiserum are comparable to those of pooled sheep serum. No foreign proteins added. It is reconstituted by adding 1 ml sterile
Source	Purified Fab of normal mouse IgG. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Delipidated, heat inactivated lyophilized stable whole serum No preservative added. Total protein and IgG concentration in the antiserum are comparable to those of pooled sheep serum. No foreign proteins added.

Application Notes

Precipitation assays. In immunoelectrophoresis use 2 μ l or equivalent against 120 μ l antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 μ l antiserum in a 3 mm diameter centre well and 2 μ l serum samples (neat and