

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Sheep anti Rat IgG Fab, Clone: [Polyclonal], Monoclonal NMB-SHARA/FAB/7S

| Article Name | Sheep anti Rat IgG Fab, Clone: [Polyclonal], Monoclonal |
|----------------------------|--|
| Biozol Catalog Number | NMB-SHARA/FAB/7S |
| Supplier Catalog Number | ShARa/Fab/7S |
| Alternative Catalog Number | NMB-SHARA/FAB/7S |
| Manufacturer | NordicMubio |
| Host | Sheep |
| Category | Antikörper |
| Product Description | The antiserum is reacting with the Fab subunit of intact IgG, IgA and IgM and other Ig classes of both light chain types, with their Fab or F(ab)2 subunits and with free light chains of kappa and lambda type as tested in immunoelectrophoresis and dou |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| Buffer | Purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1 ml sterile distilled w |
| Source | Purified Fab from normal IgG isolated from pooled rat serum. Freunds complete adjuvant is used in the first step of the immunization procedure. |
| Formula | Purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. |

Indirect immunofluorescence, ELISA, Dot blot, Immunoblotting.