

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Sheep anti Rat IgG Fab, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-SHARA/FAB/BIO

Article Name	Sheep anti Rat IgG Fab, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-SHARA/FAB/BIO
Supplier Catalog Number	ShARa/Fab/Bio
Alternative Catalog Number	NMB-SHARA/FAB/BIO
Manufacturer	NordicMubio
Host	Sheep
Category	Antikörper
Conjugation	Biotin
Product Description	The antiserum is reacting with the Fab subunit of intact IgG, IgA and IgM and other Ig classes of both light chain types, with their Fab or F(ab)2 subunits and with free light chains of kappa and lambda type as tested in immunoelectrophoresis and dou
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Biotin-coupled purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into

Source	Purified Fab from normal IgG isolated from pooled rat serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Biotin-coupled purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting.