

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Sheep anti Rat IgG Fab, conjugated with TRITC, Clone: [Polyclonal], Monoclonal NMB-SHARA/FAB/TRITC

Article Name	Sheep anti Rat IgG Fab, conjugated with TRITC, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-SHARA/FAB/TRITC
Supplier Catalog Number	ShARa/Fab/TRITC
Alternative Catalog Number	NMB-SHARA/FAB/TRITC
Manufacturer	NordicMubio
Host	Sheep
Category	Antikörper
Conjugation	TRITC
Product Description	The reactivity of the antiserum is directed to the Fab subunits of intact IgG, IgA IgM and other Ig classes of both light chain types, with their Fab or F(ab)2 subunits, and with free light chains of kappa and lambda type. In immunoelectrophoresis an
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	TRITC-coupled purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added as it may interfere with he antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spun dow
Source	Purified Fab of normal IgG isolated form rat serum. Freunds complete adjuvant is used in the first step of the immunization procedure.

Formula	TRITC-coupled purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added as it may interfere with he antibody activity.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence.