



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Swine anti Rabbit IgG (heavy and light chains), conjugated with FITC, Clone: [Polyclonal], Porcine, Monoclonal NMB-SWAR/IGG(H+L)/FITC

Article Name	Swine anti Rabbit IgG (heavy and light chains), conjugated with FITC, Clone: [Polyclonal], Porcine, Monoclonal
Biozol Catalog Number	NMB-SWAR/IGG(H+L)/FITC
Supplier Catalog Number	SwAR/IgG(H+L)/FITC
Alternative Catalog Number	NMB-SWAR/IGG(H+L)/FITC
Manufacturer	NordicMubio
Host	Porcine
Category	Antikörper
Conjugation	FITC
Product Description	The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in rabbit serum, a
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	FITC-coupled purified hyperimmune swine IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 2 ml sterile distilled water, spun do

Source	Purified normal IgG isolated from pooled rabbit serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	FITC-coupled purified hyperimmune swine IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence.