

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-GFAP Purified, Clone: [GF-02], Monoclonal EXB-11-441-C025

Article Name Anti-GFAP Purified, Clone: [GF-02], Monoclonal Biozol Catalog Number EXB-11-441-C025 Supplier Catalog Number 11-441-C025 Alternative Catalog Number EXB-11-441-C025 Manufacturer EXBIO Category Antikörper Application WB, IHC-P, IHC-Fr, ICC Species Reactivity Human, Porcine Immunogen Pellet of porcine brain cold-stable proteins after depolymerization of microtubules. GFAP (glial fibrillary acidic protein) was discovered by Bignami et al.
Supplier Catalog Number 11-441-C025 Alternative Catalog Number EXB-11-441-C025 Manufacturer EXBIO Category Antikörper Application WB, IHC-P, IHC-Fr, ICC Species Reactivity Human, Porcine Immunogen Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.
Alternative Catalog Number EXB-11-441-C025 Manufacturer EXBIO Category Antikörper Application WB, IHC-P, IHC-Fr, ICC Species Reactivity Human, Porcine Immunogen Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.
Manufacturer EXBIO Category Antikörper Application WB, IHC-P, IHC-Fr, ICC Species Reactivity Human, Porcine Immunogen Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.
Category Antikörper Application WB, IHC-P, IHC-Fr, ICC Species Reactivity Human, Porcine Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.
Application WB, IHC-P, IHC-Fr, ICC Species Reactivity Human, Porcine Immunogen Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.
Species Reactivity Human, Porcine Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.
Immunogen Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.
microtubules.
GFAP (glial fibrillary acidic protein) was discovered by Bignami et a
Product Description (1972) as a major fibrous protein of multiple sclerosis plaques. It we subsequently found to be a member of the 10 nm or intermediate filament protein family, specifically the i
Clonality Monoclonal
Concentration 1 mg/ml
Clone Designation [GF-02]
Isotype IgM
Buffer Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Storage 2°C to 8°C

Antibody Type	Monoclonal Antibody
Application Notes	Western blotting: Recommended dilution: 1-2 µg/ml.