

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Anti-Hu CD3 PerCP-Cy(TM)5.5, Clone: [SK7], PerCP/Cy5.5, Monoclonal EXB-T9-173-T100

|                            |   |
|----------------------------|---|
| Article Name               | Anti-Hu CD3 PerCP-Cy(TM)5.5, Clone: [SK7], PerCP/Cy5.5, Monoclonal  |
| Biozol Catalog Number      | EXB-T9-173-T100   |
| Supplier Catalog Number    | T9-173-T100   |
| Alternative Catalog Number | EXB-T9-173-T100   |
| Manufacturer               | EXBIO   |
| Category                   | Antikörper  |
| Application                | FC  |
| Species Reactivity         | Human, Primate  |
| Immunogen                  | Human thymocytes  |
| Conjugation                | PerCP/Cy5.5   |
| Product Description        | CD3 complex is crucial in transducing antigen-recognition signals into the cytoplasm of T cells and in regulating the cell surface expression of the TCR complex. T cell activation through the antigen receptor (TCR) involves the cytoplasmic tails of t... |
| Clonality                  | Monoclonal  |
| Clone Designation          | [SK7]   |
| Isotype                    | Mouse IgG1 kappa  |
| Buffer                     | Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide   |
| Storage                    | 2°C to 8°C  |

|                    |  |
|--------------------|--|
| Target             | CD3  |
| Antibody Type      | Monoclonal Antibody  |
| Application Dilute | Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.   |
| Application Notes  | Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.<br>Extracellular and intracellular stai |