

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Acetyltransferase Assay Buffer SCM-A01-09

Article Name	Acetyltransferase Assay Buffer
Biozol Catalog Number	SCM-A01-09
Supplier Catalog Number	A01-09
Alternative Catalog Number	SCM-A01-09-1,SCM-A01-09-5,SCM-A01-09-20
Manufacturer	SignalChem
Category	Kits/Assays
Application	Acetyltransferase Assay
Product Description	Acetyltransferase Assay Buffer...
Buffer	Frozen solution of buffer containing: 250mM Tris-HCl, pH 8.0, 0.5mM EDTA and 25% glycerol. Add to final concentration of 2mM DTT (Catalog D86-09B) prior to use. DTT Solution supplied in separate vial.
Source	5X stock solution of the buffer used for assaying histone acetyltransferases (KAT).
Form	Frozen solution of buffer containing: 250mM Tris-HCl, pH 8.0, 0.5mM EDTA and 25% glycerol. Add to final concentration of 2mM DTT (Catalog# D86-09B) prior to use. DTT Solution supplied in separate vial.
Formula	Frozen solution of buffer containing: 250mM Tris-HCl, pH 8.0, 0.5mM EDTA and 25% glycerol. Add to final concentration of 2mM DTT (Catalog D86-09B) prior to use. DTT Solution supplied in separate vial.