

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Acetyltransferase Dilution Buffer SCM-A21-09

Article Name	Acetyltransferase Dilution Buffer
Biozol Catalog Number	SCM-A21-09
Supplier Catalog Number	A21-09
Alternative Catalog Number	SCM-A21-09-1,SCM-A21-09-5,SCM-A21-09-20
Manufacturer	SignalChem
Category	Sonstiges
Application	Acetyltransferase Assay
Product Description	Acetyltransferase Dilution Buffer...
Buffer	Frozen solution of buffer containing: 50mM Tris-HCl, pH 8.0, 0.1mM EDTA and 5% glycerol. Add to final concentration of 0.4mM DTT (Catalog D86-09B) prior to use. DTT Solution supplied in separate vial.
Source	Solution of buffer used for diluting enzyme histone acetyltransferases (KAT).
Form	Frozen solution of buffer containing: 50mM Tris-HCl, pH 8.0, 0.1mM EDTA and 5% glycerol. Add to final concentration of 0.4mM DTT (Catalog# D86-09B) prior to use. DTT Solution supplied in separate vial.
Formula	Frozen solution of buffer containing: 50mM Tris-HCl, pH 8.0, 0.1mM EDTA and 5% glycerol. Add to final concentration of 0.4mM DTT (Catalog D86-09B) prior to use. DTT Solution supplied in separate vial.