



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

AMP Solution SCM-A46-09

Article Name	AMP Solution
Biozol Catalog Number	SCM-A46-09
Supplier Catalog Number	A46-09
Alternative Catalog Number	SCM-A46-09-500
Manufacturer	SignalChem
Category	Sonstiges
Application	Kinase Assay
Product Description	AMP Solution
Buffer	Frozen solution of buffer containing: 25mM MOPS pH 7.2, and 0.5mM AMP.
Source	0.5mM AMP solution used for AMPK family kinase assays.
Form	Frozen solution of buffer containing: 25mM MOPS pH 7.2, and 0.5mM AMP.
Formula	Frozen solution of buffer containing: 25mM MOPS pH 7.2, and 0.5mM AMP.