

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### ATP Stock Solution SCM-A50-09

Article Name	ATP Stock Solution
Biozol Catalog Number	SCM-A50-09
Supplier Catalog Number	A50-09
Alternative Catalog Number	SCM-A50-09-200
Manufacturer	SignalChem
Category	Sonstiges
Application	Kinase Assay
Product Description	ATP Stock Solution...
Buffer	Frozen solution containing: 10 mM ATP, 25 mM MOPS, pH 7.2, 12.5mM B-glycerol-phosphate, 25mM MgCl <sub>2</sub> , 5mM EGTA, 2mM EDTA, 0.25mM DTT.
Source	10 mM ATP Stock Solution used for preparing hot ATP assay cocktail.
Form	Frozen solution containing: 10 mM ATP, 25 mM MOPS, pH 7.2, 12.5mM B-glycerol-phosphate, 25mM MgCl <sub>2</sub> , 5mM EGTA, 2mM EDTA, 0.25mM DTT.
Formula	Frozen solution containing: 10 mM ATP, 25 mM MOPS, pH 7.2, 12.5mM B-glycerol-phosphate, 25mM MgCl <sub>2</sub> , 5mM EGTA, 2mM EDTA, 0.25mM DTT.