

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Kinase Dilution Buffer X SCM-K20-09

Article Name	Kinase Dilution Buffer X
Biozol Catalog Number	SCM-K20-09
Supplier Catalog Number	K20-09
Alternative Catalog Number	SCM-K20-09-1,SCM-K20-09-5,SCM-K20-09-20
Manufacturer	SignalChem
Category	Proteine/Peptide
Application	Kinase Assay
Product Description	Kinase Dilution Buffer X...
Buffer	Frozen solution of buffer containing: 40mM Tris-HCl, pH7.4, 20mM MgCl ₂ , 2.5mM MnCl ₂ , and 0.1mg/ml BSA.
Source	The working solution used for diluting kinase proteins.
Form	Frozen solution of buffer containing: 40mM Tris-HCl, pH7.4, 20mM MgCl ₂ , 2.5mM MnCl ₂ , and 0.1mg/ml BSA.
Formula	Frozen solution of buffer containing: 40mM Tris-HCl, pH7.4, 20mM MgCl ₂ , 2.5mM MnCl ₂ , and 0.1mg/ml BSA.