

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Kinase Dilution Buffer II SCM-K22-09

Article Name	Kinase Dilution Buffer II
Biozol Catalog Number	SCM-K22-09
Supplier Catalog Number	K22-09
Alternative Catalog Number	SCM-K22-09-1,SCM-K22-09-5,SCM-K22-09-20
Manufacturer	SignalChem
Category	Proteine/Peptide
Application	Kinase Assay
Product Description	Kinase Dilution Buffer II...
Buffer	Frozen solution of buffer containing: 5 mM MOPS, pH7.2, 2.5mM B-glycerol-phosphate, 4mM MgCl <sub>2</sub> , 2.5mM MnCl <sub>2</sub> , 1mM EGTA, 0.4mM EDTA, 0.05mM DTT.
Source	Solution of buffer used for diluting kinase protein.
Form	Frozen solution of buffer containing: 5 mM MOPS, pH7.2, 2.5mM B-glycerol-phosphate, 4mM MgCl <sub>2</sub> , 2.5mM MnCl <sub>2</sub> , 1mM EGTA, 0.4mM EDTA, 0.05mM DTT.
Formula	Frozen solution of buffer containing: 5 mM MOPS, pH7.2, 2.5mM B-glycerol-phosphate, 4mM MgCl <sub>2</sub> , 2.5mM MnCl <sub>2</sub> , 1mM EGTA, 0.4mM EDTA, 0.05mM DTT.