



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Lipid Dilution Buffer II SCM-L22-09

Article Name	Lipid Dilution Buffer II
Biozol Catalog Number	SCM-L22-09
Supplier Catalog Number	L22-09
Alternative Catalog Number	SCM-L22-09-1,SCM-L22-09-5,SCM-L22-09-20
Manufacturer	SignalChem
Category	Sonstiges
Application	Kinase Assay
Product Description	Lipid Dilution Buffer II
Buffer	Frozen solution of buffer containing: 25mM HEPES, pH 7.5, 0.5mM EGTA, and 0.5% Triton X-100.
Source	The working solution used for diluting lipid substrates.
Form	Frozen solution of buffer containing: 25mM HEPES, pH 7.5, 0.5mM EGTA, and 0.5% Triton X-100.
Formula	Frozen solution of buffer containing: 25mM HEPES, pH 7.5, 0.5mM EGTA, and 0.5% Triton X-100.