



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## EnzyChrom(TM) Glycogen Assay Kit BAS-E2GN-100

Article Name	EnzyChrom(TM) Glycogen Assay Kit
Biozol Catalog Number	BAS-E2GN-100
Supplier Catalog Number	E2GN-100
Alternative Catalog Number	BAS-E2GN-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of glycogen and evaluation of drug effects on glycogen metabolism. Key Features: Use as little as 10 $\mu$ L samples. Linear detection range: 2 to 200 $\mu$ g/mL glycogen for colorimetric assays and 0.2 to 20 $\mu$ g/mL for fluorimetr