



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## EnzyChrom(TM) NAD/NADH Assay Kit BAS-E2ND-100

| Article Name               | EnzyChrom(TM) NAD/NADH Assay Kit  |
|----------------------------|---|
| Biozol Catalog Number      | BAS-E2ND-100  |
| Supplier Catalog Number    | E2ND-100  |
| Alternative Catalog Number | BAS-E2ND-100  |
| Manufacturer               | BioAssay Systems  |
| Category                   | Kits/Assays   |
| Product Description        | For sensitive determination of NAD and NADH and evaluation of drug effects on NAD/NADH metabolism. Key Features: Sensitive and accurate. Detection limit of 0.05 $\mu\text{M}$ and linearity up to 10 $\mu\text{M}$ NAD+/NADH in 96-well plate assay. Convenient. The procedure |